

MidAtlantic Diagnostics Announces the

IVF-1 CHAMBER

Mobile • Self-Contained • IVF Workstation • One-Year Full Warranty • Refurbished Units Available



STANDARD FEATURES

User-Defined Temperature Setting
Integrated microscope (new or existing)
External Fiber optic light source
CO2 Controller and Regulator

CUSTOM SOLUTIONS & OPTIONS

Unique and standard microscopes
Integrated Battery Backup
Integrated Video setups available
Integrated Humidification System

MidAtlantic
DIAGNOSTICS, INC.

438 N. Elmwood Road Marlton, NJ 08053 • Phone: 856-762-2000 • 800-648-1151 • Fax: 856-762-2009

www.midatlanticdiagnostics.com

IVF-1 Chamber (Hoffman) Accessories and Options

Battery Back Up

Description

Integrated system for 60 minutes uninterrupted power supply in event of power outage. System will automatically switch to battery power and automatically switch back to primary AC power when available.

Part Number

CP-IVF1-BAT1200

Integrated system for 30 minutes uninterrupted power supply in event of power outage. System will automatically switch to battery power and automatically switch back to primary AC power when available.

CP-IVF1-BAT900

Elevation Inserts

Description

Mount between cabinet body and base to elevate the working height of the IVF-1. Available in five sizes:

Two-inch insert, Elevates unit two inches.

Three-inch insert, Elevates unit three inches.

Four-inch insert, Elevates unit four inches.

Five-inch insert, Elevates unit five inches.

Six-inch insert, Elevates unit six inches.

Part Number

CP-IVF1-Insert2

CP-IVF1-Insert3

CP-IVF1-Insert4

CP-IVF1-Insert5

CP-IVF1-Insert6

Internal Shelf System

Description

Five (5) tiers of shelves for temporary storage of culture dishes.

Shelf mounts in either right or left corners of IVF-1.

Easily removed for cleaning.

Part Number

CP-IVF1-Shelf

Spacer and Insulation Block

Description

Adds an insulation space between the work surface around the scope and the heat source under scope. Recommended if material is left on work surface for extended period.

1-inch Spacer

1.25-inch Spacer

1.50-inch Spacer

2-inch Spacer

2.25-inch Spacer

2.50-inch Spacer

Other custom sizes available upon request

Part Number

CP-IVF1-Spacer1

CP-IVF1-Spacer125

CP-IVF1-Spacer150

CP-IVF1-Spacer2

CP-IVF1-Spacer225

CP-IVF1-Spacer250

Preventative Maintenance Kits

Description

Replacement parts kit for IVF-1 including gasket, air filter, tubing, and Forma filter and sleeves.

For IVF-1 with OVAL access ports

For IVF-1 with ROUND access ports

Part Number

CP-IVF1-PMOval

CP-IVF1-PMRound

