





Assisted Hatching


Hatching Overall Length




MAH-SM



MAH-LG



Hatching Angle Options



15°, 20°, 25°, 30°, 35°, 40°, 45°

Hatching Side View



Straight

Assisted Hatching (Chemical) micropipets are used to create a hole in the zona pellucida using acid Tyrode's solution. These micropipets are made with a blunt opening, slightly polished, and are offered in two sizes. Other diameters may be ordered as customized Assisted Hatching micropipets by specifying the desired size of the opening.

Hatching Micropipets (10 per box)

Reference No.	Description	I.D.
MAH-SM-0	Straight Hatching	8-10 μ m
MAH-SM-15	15° angled Hatching	8-10 μ m
MAH-SM-20	20° angled Hatching	8-10 μ m
MAH-SM-25	25° angled Hatching	8-10 μ m
MAH-SM-30	30° angled Hatching	8-10 μ m
MAH-SM-35	35° angled Hatching	8-10 μ m
MAH-SM-40	40° angled Hatching	8-10 μ m
MAH-SM-45	45° angled Hatching	8-10 μ m
Reference No.	Description	I.D.
MAH-LG-0	Straight Hatching	10-12 μ m
MAH-LG-15	15° angled Hatching	10-12 μ m
MAH-LG-20	20° angled Hatching	10-12 μ m
MAH-LG-25	25° angled Hatching	10-12 μ m
MAH-LG-30	30° angled Hatching	10-12 μ m
MAH-LG-35	35° angled Hatching	10-12 μ m
MAH-LG-40	40° angled Hatching	10-12 μ m
MAH-LG-45	45° angled Hatching	10-12 μ m