





# Holding Micropipets

Overall Length




Holding Diameters


MPH-SM



MPH-MED



Angle Options



15°, 20°, 25°, 30°, 35°, 40°, 45°

Holding Side View



Straight

**Holding micropipets** are used to hold and immobilize the oocyte or blastocyst. Standard Humagen Holding micropipets have a blunt, polished opening. Customized Holding micropipets may also be ordered by specifying an exact size (O.D. or I.D.) rather than the indicated ranges.

## Small Holding Micropipets (10 per box)



Reference No.	Description	I.D.	O.D.
MPH-SM-0	Straight Small Holding	15-20µm	65-95µm
MPH-SM-15	15° angled Small Holding	15-20µm	65-95µm
MPH-SM-20	20° angled Small Holding	15-20µm	65-95µm
MPH-SM-25	25° angled Small Holding	15-20µm	65-95µm
MPH-SM-30	30° angled Small Holding	15-20µm	65-95µm
MPH-SM-35	35° angled Small Holding	15-20µm	65-95µm
MPH-SM-40	40° angled Small Holding	15-20µm	65-95µm
MPH-SM-45	45° angled Small Holding	15-20µm	65-95µm

## Medium Holding Micropipets (10 per box)

Reference No.	Description	I.D.	O.D.
MPH-MED-0	Straight Medium Holding	15-20µm	95-120µm
MPH-MED-15	15° angled Medium Holding	15-20µm	95-120µm
MPH-MED-20	20° angled Medium Holding	15-20µm	95-120µm
MPH-MED-25	25° angled Medium Holding	15-20µm	95-120µm
MPH-MED-30	30° angled Medium Holding	15-20µm	95-120µm
MPH-MED-35	35° angled Medium Holding	15-20µm	95-120µm
MPH-MED-40	40° angled Medium Holding	15-20µm	95-120µm
MPH-MED-45	45° angled Medium Holding	15-20µm	95-120µm

### Large Holding Micropipets (10 per box)


Reference No.	Description	I.D.	O.D.
MPH-LG-0	Straight Large Holding	25-30 $\mu$ m	120-150 $\mu$ m
MPH-LG-15	15° angled Large Holding	25-30 $\mu$ m	120-150 $\mu$ m
MPH-LG-20	20° angled Large Holding	25-30 $\mu$ m	120-150 $\mu$ m
MPH-LG-25	25° angled Large Holding	25-30 $\mu$ m	120-150 $\mu$ m
MPH-LG-30	30° angled Large Holding	25-30 $\mu$ m	120-150 $\mu$ m
MPH-LG-35	35° angled Large Holding	25-30 $\mu$ m	120-150 $\mu$ m



### Extra Large Holding Micropipets (10 per box)

Reference No.	Description	I.D.	O.D.
MPH-XLG-0	Straight Extra Large Holding	25-30 $\mu$ m	150-180 $\mu$ m
MPH-XLG-15	15° angled X-Large Holding	25-30 $\mu$ m	150-180 $\mu$ m
MPH-XLG-20	20° angled X-Large Holding	25-30 $\mu$ m	150-180 $\mu$ m
MPH-XLG-25	25° angled X-Large Holding	25-30 $\mu$ m	150-180 $\mu$ m
MPH-XLG-30	30° angled X-Large Holding	25-30 $\mu$ m	150-180 $\mu$ m
MPH-XLG-35	35° angled X-Large Holding	25-30 $\mu$ m	150-180 $\mu$ m

### Angle Options



15°, 20°, 25°, 30°, 35°