



Standard ICSI

ICSI Angle Options



15°, 20°, 25°, 30°, 35°, 40°, 45°

Overall Length




ICSI Side View



Side View

With Spike



MIC-50






MIC-SI



MIC-35

Front View



ICSI Micropipets

ICSI (Intracytoplasmic Sperm Injection) micropipets are used to aspirate and inject the sperm cell directly into the oocyte. **Larger ICSI micropipets** are used for aspiration and injecting immature sperm. The 8-9µm size is a standard Spermatid micropipet. Other diameters are available by ordering a customized ICSI micropipet. Available options without customizing: change overall length, change distance of bend from tip, no tip protectors, no spike on bevel.

ICSI Micropipets with Spike (10 per box)

Reference No.	Description	Bevel Length	I.D.
MIC-50-0	Straight ICSI	8-10µm	5-6µm
MIC-50-15	15° angled ICSI	8-10µm	5-6µm
MIC-50-20	20° angled ICSI	8-10µm	5-6µm
MIC-50-25	25° angled ICSI	8-10µm	5-6µm
MIC-50-30	30° angled ICSI	8-10µm	5-6µm
MIC-50-35	35° angled ICSI	8-10µm	5-6µm
MIC-50-40	40° angled ICSI	8-10µm	5-6µm
MIC-50-45	45° angled ICSI	8-10µm	5-6µm


Reference No.	Description	Bevel Length	I.D.
MIC-SI-0	Straight ICSI	10-11µm	5-6µm
MIC-SI-15	15° angled ICSI	10-11µm	5-6µm
MIC-SI-20	20° angled ICSI	10-11µm	5-6µm
MIC-SI-25	25° angled ICSI	10-11µm	5-6µm
MIC-SI-30	30° angled ICSI	10-11µm	5-6µm
MIC-SI-35	35° angled ICSI	10-11µm	5-6µm
MIC-SI-40	40° angled ICSI	10-11µm	5-6µm
MIC-SI-45	45° angled ICSI	10-11µm	5-6µm

Reference No.	Description	Bevel Length	I.D.
MIC-35-0	Straight ICSI	11-12µm	5-6µm
MIC-35-15	15° angled ICSI	11-12µm	5-6µm
MIC-35-20	20° angled ICSI	11-12µm	5-6µm
MIC-35-25	25° angled ICSI	11-12µm	5-6µm
MIC-35-30	30° angled ICSI	11-12µm	5-6µm
MIC-35-35	35° angled ICSI	11-12µm	5-6µm
MIC-35-40	40° angled ICSI	11-12µm	5-6µm
MIC-35-45	45° angled ICSI	11-12µm	5-6µm

ICSI Micropipets with Spike (10 per box)

Reference No.	Description	Bevel Length	I.D.
MIC-6-0	Straight ICSI	9-10 μ m	4-5 μ m
MIC-6-15	15° angled ICSI	9-10 μ m	4-5 μ m
MIC-6-20	20° angled ICSI	9-10 μ m	4-5 μ m
MIC-6-25	25° angled ICSI	9-10 μ m	4-5 μ m
MIC-6-30	30° angled ICSI	9-10 μ m	4-5 μ m
MIC-6-35	35° angled ICSI	9-10 μ m	4-5 μ m
MIC-6-40	40° angled ICSI	9-10 μ m	4-5 μ m
MIC-6-45	45° angled ICSI	9-10 μ m	4-5 μ m


MIC-6




Larger ICSI Micropipets (10 per box)

Reference No.	Description	Bevel Length	I.D.
MIC-8-0	Straight ICSI	8-10 μ m	7-8 μ m
MIC-8-15	15° angled ICSI	8-10 μ m	7-8 μ m
MIC-8-20	20° angled ICSI	8-10 μ m	7-8 μ m
MIC-8-25	25° angled ICSI	8-10 μ m	7-8 μ m
MIC-8-30	30° angled ICSI	8-10 μ m	7-8 μ m
MIC-8-35	35° angled ICSI	8-10 μ m	7-8 μ m
MIC-8-40	40° angled ICSI	8-10 μ m	7-8 μ m
MIC-8-45	45° angled ICSI	8-10 μ m	7-8 μ m

MIC-8



MIC-9



Standard Spermatid ICSI Micropipets (10 per box)

Reference No.	Description	Bevel Length	I.D.
MIC-9-0	Straight ICSI	9-11 μ m	8-9 μ m
MIC-9-15	15° angled ICSI	9-11 μ m	8-9 μ m
MIC-9-20	20° angled ICSI	9-11 μ m	8-9 μ m
MIC-9-25	25° angled ICSI	9-11 μ m	8-9 μ m
MIC-9-30	30° angled ICSI	9-11 μ m	8-9 μ m
MIC-9-35	35° angled ICSI	9-11 μ m	8-9 μ m
MIC-9-40	40° angled ICSI	9-11 μ m	8-9 μ m
MIC-9-45	45° angled ICSI	9-11 μ m	8-9 μ m