




PZD Micropipets

PZD (Partial Zona Dissection) micropipets are used in blastomere biopsy and mechanical hatching techniques by physically piercing the zona pellucida. Humagen PZD micropipets are pulled with a long thin taper to a sharp point.

PZD Micropipets (10 per box)

Reference No.	Description
PZD-0	Straight PZD micropipet
PZD-15	15° angled PZD micropipet
PZD-20	20° angled PZD micropipet
PZD-25	25° angled PZD micropipet
PZD-30	30° angled PZD micropipet
PZD-35	35° angled PZD micropipet
PZD-40	40° angled PZD micropipet
PZD-45	45° angled PZD micropipet


PZD Overall Length



Straight PZD



PZD Angle Options



15°, 20°, 25°, 30°, 35°, 40°, 45°