





# Polar Body Biopsy

**Polar Body Biopsy micropipets** (for PGD) are made with a 13-15 $\mu$ m I.D. opening and are offered flat and polished, beveled and polished or beveled with a spike. Other diameters are available by ordering customized biopsy micropipets and specifying the desired diameter.


Polar Body Overall Length




MPB-FP




MPB-BP



MPB-BS



Polar Body Angle Options



15°, 20°, 25°, 30°, 35°, 40°, 45°

## Polar Body Biopsy Micropipets (10 per box)

Reference No.	Description	I.D.	Spike
MPB-FP-0	Straight, flat and polished	13-15 $\mu$ m	No
MPB-FP-15	15° angled, flat and polished	13-15 $\mu$ m	No
MPB-FP-20	20° angled, flat and polished	13-15 $\mu$ m	No
MPB-FP-25	25° angled, flat and polished	13-15 $\mu$ m	No
MPB-FP-30	30° angled, flat and polished	13-15 $\mu$ m	No
MPB-FP-35	35° angled, flat and polished	13-15 $\mu$ m	No
MPB-FP-40	40° angled, flat and polished	13-15 $\mu$ m	No
MPB-FP-45	45° angled, flat and polished	13-15 $\mu$ m	No
Reference No.	Description	I.D.	Spike
MPB-BP-0	Straight, beveled and polished	13-15 $\mu$ m	No
MPB-BP-15	15° angled, beveled and polished	13-15 $\mu$ m	No
MPB-BP-20	20° angled, beveled and polished	13-15 $\mu$ m	No
MPB-BP-25	25° angled, beveled and polished	13-15 $\mu$ m	No
MPB-BP-30	30° angled, beveled and polished	13-15 $\mu$ m	No
MPB-BP-35	35° angled, beveled and polished	13-15 $\mu$ m	No
MPB-BP-40	40° angled, beveled and polished	13-15 $\mu$ m	No
MPB-BP-45	45° angled, beveled and polished	13-15 $\mu$ m	No
Reference No.	Description	I.D.	Spike
MPB-BS-0	Straight, beveled with spike	13-15 $\mu$ m	Yes
MPB-BS-15	15° angled, beveled with spike	13-15 $\mu$ m	Yes
MPB-BS-20	20° angled, beveled with spike	13-15 $\mu$ m	Yes
MPB-BS-25	25° angled, beveled with spike	13-15 $\mu$ m	Yes
MPB-BS-30	30° angled, beveled with spike	13-15 $\mu$ m	Yes
MPB-BS-35	35° angled, beveled with spike	13-15 $\mu$ m	Yes
MPB-BS-40	40° angled, beveled with spike	13-15 $\mu$ m	Yes
MPB-BS-45	45° angled, beveled with spike	13-15 $\mu$ m	Yes