




Denuding Micropipets

Denuding Overall Length



- 120 μm I.D.
- 130 μm I.D.
- 135 μm I.D.
- 140 μm I.D.
- 150 μm I.D.
- 160 μm I.D.
- 170 μm I.D.
- 175 μm I.D.
- 190 μm I.D.

Denuding micropipets are used to remove the cumulus cells from the oocyte. These micropipets are made from 1.5 mm O.D. glass, have a blunt opening and are sold straight. They are available in a variety of sizes. Denudation handle sold separately. Silicon bulbs included with pipets.

Denuding Micropipets (10 per box)

Reference No.	Description	I.D.
MDP-120-0	Straight Denuding Micropipet, with blunt opening	117-123 μm
MDP-130-0	Straight Denuding Micropipet, with blunt opening	127-133 μm
MDP-135-0	Straight Denuding Micropipet, with blunt opening	132-138 μm
MDP-140-0	Straight Denuding Micropipet, with blunt opening	137-143 μm
MDP-150-0	Straight Denuding Micropipet, with blunt opening	147-153 μm
MDP-160-0	Straight Denuding Micropipet, with blunt opening	157-163 μm
MDP-170-0	Straight Denuding Micropipet, with blunt opening	167-173 μm
MDP-175-0	Straight Denuding Micropipet, with blunt opening	172-178 μm
MDP-190-0	Straight Denuding Micropipet, with blunt opening	187-193 μm

Denuding Handle

Reference No.	Description
HANDLE-BL	Black Denuding Handle

