

# XT Series



## Features

- New Rugged Construction – ribbed high strength aluminum body, magniformed necktube design, and more durable paint
- Versatile storage system - convenient canister index location ring and internal spider.
- Superior vacuum performance with super insulation provides maximum holding times
- Security - Accessory low-level alarm is available for more valuable samples
- Mobility - Roller bases are available

## XT Series

The Extended Time Series (XT) of cryogenic refrigerators are designed for long-term storage of a variety of materials at cryogenic temperatures. This series also offers a low profile model (XTL) for use in more confined workspaces.

Models	XTL3	XTL8	XT10	XT20	XT34	
Static Holding Time days (1)	27	80	100	230	340	
Working Time days (2)	17	50	62	140	212	
Evaporation Rate (1) liters/days	0.11	0.10	0.10	0.09	0.10	
Liquid Nitrogen Capacity - liters	3	8	10	20.7	34	
Weight Empty	- lbs.	7.2	19.6	16.5	26	34.75
	- kg.	3.3	8.9	7.5	11.8	15.8
Weight Full (3)	- lbs.	12.5	33.8	34.3	62.8	95.3
	- kg.	5.7	15.4	15.6	28.6	43.3
Neck Diameter	- in.	2.0	2.0	2.0	2.0	2.0
	- mm.	51	51	51	51	51
Overall Height	- in.	17.2	19	23.8	25.8	26.3
	- mm.	437	483	597	655	668
Overall Diameter	- in.	7.6	15.6	11.4	15.6	18.8
	- mm.	193	396	290	396	478
Number of Canisters	6	6	6	6	6	
Canister Dimensions (4) - in.	1.5x5	1.5x5	1.5x11	1.5x11	1.5x11	
	- mm.	38x127	38.127	38x279	38x279	38x279
Number of 1.2 ml & 2.0 ml vials (5/cane)	NA	NA	150	150	150	
Number of 1.2 ml & 2.0 ml vials (6/cane)	NA	NA	180	180	180	
Number of 1/2 cc straws (10/cane)	NA	NA	540	540	540	
Number of 1.2 cc straws - Bulk (1 level)	750	750	750	750	750	
Number of 1/2 cc straws - Bulk (2 levels)	NA	NA	1500	1500	1500	
Roller Base	N/A	R018-8C00	N/A	R018-8C00	R033-8C00	
Low-Level Alarm	N/A	N/A	N/A	R033-8C15	R033-8C15	

(1) Evaporation rate and static holding time are nominal. Actual rate may be affected by the nature of the contents, atmospheric conditions, container history, and manufacturing tolerances.

(2) Work time is an arbitrary, reference-only value to estimate container performance under the actual operating conditions.

(3) Without canisters or racks.

(4) Canisters also available in 5.0 in. (127 mm) height for 10 XT, 20XT and 34XT.