




Vibration Isolation Platform VIP-3000 Series



Vistek Isolation The Elegant Vibration Solution

Vistek Inc.
4405 E. Baseline Rd Suite 120
Phoenix, Arizona 85042
Phone (602) 426-1211
Fax (602) 426-1611
www.vistekinc.com

A DIVISION OF



TRUE GRAVITY ENTERPRISES

Vibration Isolation Platform

VIP-3000 Series



Low profile, countertop vibration isolation platform for research microscopes and imaging systems. The VIP is comprised of a solid platform and three Vistek* Series 300 isolation bearings.

U.S. Pat. Nos. 6,517,060/6,520,283 & US and Intl. Pats. Pend.

Simple design...exceptional performance!

CONVENIENCE

- Easy installation
- Maneuverability
- Light weight

SIZE

- Compact
- Portable
- Low profile

COST

- Low purchase price
- No upkeep expenses
- Low shipping cost

MAINTENANCE FREE

- Passive mechanical
- No air hoses
- No electrical hardware

Contact your local microscope dealer for a demonstration and pricing information.
Visit our web site for a complete Vistek product line.



www.vistekinc.com