

# *MidAtlantic Diagnostics*

## Custom Products **IVF WORKSTATION**

The ORIGINAL Workstation for Your IVF Laboratory



Shown is the Six Foot IVF Workstation (Dual Scope Model)  
Also available are the 3, 4, 5 and 6 Foot IVF Workstations

### STANDARD FEATURES

Vertical Laminar Flow  
Integrated Warming Plate  
Water Circulator and Heater  
Stereo Scope Installation Included  
Quietest Performance in Industry  
HEPA Air Filtration  
Available in 3, 4, 5 & 6 Foot Models  
Dual Electrical Receptacle (110-Volt) Mounted  
on Back Wall of Work Area

### CUSTOM SOLUTIONS & OPTIONS

Compatible with any Microscope  
Custom Heating Patterns  
Complete Video Systems Available/Installed  
CO2 Manifold/Humidifier Systems  
- Single Flowmeter  
- Double Flowmeter  
- Triple Flowmeter  
Warming Blocks for Tubes, Dishes  
Mini-Incubator Chambers (Glass)

**MidAtlantic**  
DIAGNOSTICS, INC.

# Custom Products IVF WORKSTATION

The ORIGINAL Workstation for Your IVF Laboratory

## Standard Models

<u>Ordering Number</u>	<u>Description</u>	<u>Dimensions</u>
WKST-24X36	3-Foot IVF Workstation	38x78.7x28.2 inches <u>Work Area</u> 35.5x27.8x24 inches
<u>Ordering Number</u> WKST-24X48	<u>Description</u> 4-Foot IVF Workstation	<u>Dimensions</u> 49.8x78.7x28.2 inches <u>Work Area</u> 47.5x27.8x24 inches
<u>Ordering Number</u> WKST-24X60	<u>Description</u> 5-Foot IVF Workstation	<u>Dimensions</u> 61.6x78.7x28.2 inches <u>Work Area</u> 59x27.8x24 inches
<u>Ordering Number</u> WKST-24X72	<u>Description</u> 6-Foot IVF Workstation	<u>Dimensions</u> 73.48x78.7x28.2 inches <u>Work Area</u> 70.8x27.8x24 inches
<u>Ordering Number</u> WKST-24X72D	<u>Description</u> 6-Foot Dual IVF Workstation	<u>Dimensions</u> 73.48x78.7x28.2 inches <u>Work Area</u> 70.8x27.8x24 inches

## Each Custom Products IVF WORKSTATION includes:

- Stand
- Stainless Steel Heated Plate
- Mounting Assembly for One Microscope
- Water Circulator and Heater
- Hoses, Tubing and Clamps
- One-Year Warranty, All Parts & Labor

## Optional Items:

- Manifold and Humidifiers for CO2 Gas Mixture  
Pre-mounted in rear of work area, complete with jar and stopper cap.  
Placement can be customized by customer.
- WKST-flow1           Single Flowmeter Port, Single-Jar
- WKST-flow2           Double Flowmeter Port, Double-Jar
- WKST-flow3           Triple Flowmeter Port, Triple-Jar
- UV Light with Integrated Timer System

